

SPECIFICATION

PRIMUS GENERAL PURPOSE – Medium Steam Sterilizers

These sterilizers are designed for laboratory applications including research laboratories, bio-containment and animal care facilities. PRIMUS sterilizers are designed and manufactured in the USA according to Quality Management Systems which are in compliance with ISO 9001:2000, ISO 13485:2003, CMDCAS and FDA's Good Manufacturing Practice (GMP) for Medical Devices: General Regulation (21CFR Part 820). Medium size sterilizers are one segment of 19 standard sizes manufactured by PRIMUS. We can also custom manufacture a sterilizer to meet your specifications.



26" x 36" x 48"

DATA SHEET

STERILIZER DESIGN

The pressure vessel is the heart of any sterilizer. PRIMUS offers simple, straight forward, design. All PRIMUS pressure vessels are constructed of solid stainless steel and are fabricated by PRIMUS in our quality controlled ASME facility. The interior chamber surface is polished to a mirror finish of $<10\mu$ in. Ra. The brilliant Pri-Mirror[®] reflective chamber finish sets the highest standard for cleanliness and offers easy spill clean-up and resistance to staining and degradation.

PRIMUS rectangular chambers eliminate the wasted space and high utility costs common to old-style cylindrical or elliptical chambers and jackets. This provides savings in operating cost and service.

All medium sizes of PRIMUS sterilizers have pneumatically powered horizontal operating doors, designed to be efficient, reliable and inherently safe. To close the door, press and hold the door

button, if button is released before door is fully closed the door opens. An added safety feature, the doors stop automatically if an obstacle is encountered.

STERILIZER CYCLES

All PRIMUS steam sterilizers include Gravity, Vacuum and Liquid cycle settings. A Test/Vac cycle is provided for utilizing Bowie-Dick[™] tests.

General Purpose sterilizers provide for steam sterilization and decontamination of laboratory research and animal care products. Applications include wrapped and unwrapped hard goods, animal feed and bedding, cages and textiles. Liquids in self-venting and open containers are sterilized using Liquid cycle settings.

The time-tested and reliable PRIMUS PSS500 Microcomputer Control is an industrial microprocessor providing accuracy and automation for all customer requirements. This simple, versatile, control has become the first choice for both sterilizer users and service technicians.



PRIMUS Sterilizer Co., Inc. • 117 South 25th Street • Omaha, NE 68131

PURCHASING SPECIFICATION DATA SHEET

SPECIFICATION DATA SHEET

GENERAL (Options italicized)

All models include Vacuum, Gravity and Liquids cycles. *A Low Temperature Flowing Steam option provides flowing steam at low temperatures for Research laboratories for specific applications.*

VESSEL MATERIAL AND CONSTRUCTION

The sterilizer will have a double-wall design providing for a fully jacketed pressure vessel. The chamber cross section will be rectangular with width and height dimensions to include:

E	26"x36"x39"	G.1	32"x36"x48"
F	26"x36"x48"	H	36"x42"x60"
G	26"x36"x60"	I	36"x42"x84"

The vessel material will be non-laminated solid stainless steel plate/sheet. The chamber, head-ring and door material will be 316L stainless steel with an interior Pri-Mirror[®] finish of <10µ in. Ra.

DESIGN and MANUFACTURING APPROVALS

The vessel will be designed and constructed to ASME code Section VIII Division I and pressure rated for 45 PSIG and full vacuum. ASME Code Stamp and U-1 form will bear the name of PRIMUS as the specified US manufacturer.

The sterilizer will be manufactured according to Quality Management Systems which are in compliance with ISO 9001:2000, ISO 13485:2003, CMDCAS and will be UL listed and in conformance with CSA requirements.

DOOR CONSTRUCTION

A pneumatically operated horizontally sliding door will be operable from the touchpad push button. To close door, depress/hold door button. If door button is released while door is closing, door reverses direction and opens. In the open position, the insulated door will be secured behind a panel preventing the operator from contact with hot surfaces. An interlock will prevent cycle start unless the door is fully closed and secured.

GENERAL CONSTRUCTION

The vessel will be insulated with 1" semi-rigid high temperature fiberglass board/blanket insulation overlaid with formed aluminum paneling and mounted in a structural steel frame. The frame will be enamel coated and fitted with adjustable legs on self-centering floor pads. The fascia and side panels will be 16-gauge, type 304 stainless-steel, removable for easy service access.

WATER CONSERVATION

To conserve water and assure effluent exhausted to drain is 140° F or below, PRIMUS provides quench water on demand for exhaust effluent.

EQUIPMENT WARRANTY

Sterilizer pressure vessels manufactured by PRIMUS are warranted against defects in workmanship and materials under normal use and operation for fifteen years where the sterilizer is continually maintained under PRIMUS service contract.

SERVICE and EQUIPMENT ACCESS

Standard service access, when facing the unit, will be from left side and top. Wiring will be laid side-by-side and mechanically secured flat against the metal insulation cover. All wiring will be clearly labeled or readily visible for visual tracing. Piping components will be threaded rigid brass and flared copper fittings, positioned with sufficient space for removal and replacement without disassembly of the entire piping assembly. Wiring and piping components will be non-proprietary, industrial grade, available through Authorized Service Agencies, local supply house, or direct from PRIMUS.

STEAM SOURCE

Steam will be from an in-house steam source in the quantity and quality specified. *Optional electric steam boiler, clean steam boiler or steam-to-steam generator may be specified. Contact PRIMUS for specific boiler information.*

CONTROLS

A PSS500 microcomputer will offer a selection of eight (8) separate programs for the sterilization of wrapped goods, hard-goods, liquids and test cycles. An optional automatic Vacuum Leak Test cycle is provided to ensure the seal integrity of the vessel and piping.

Cycle parameters may be set by the user and are adaptable to a wide range of products. Alarms are annunciated, displayed and printed. Cycle process parameters will be displayed and recorded throughout the cycle. A side-mounted control panel and a thermal printer dot-matrix technology, and 32 characters per line printing is provided. A copy of the sterilization cycle report will be available at the end of the run.





PRIMUS PSS500 and PRI-MATIC[®] CONTROL PANELS WITH PRINTER



PRIMUS General Purpose - Medium Sized Autoclaves

SINGLE DOOR



VOLUME/DIMENSIONS CHART

*Refer to General Arrangement (GA) drawings for details and final connection point to utility services (S-Steam, W-Water, D-Drain, E-Electrical, A-Air).

MODEL	E	F	G	G1	H	I
Chamber Size (w x h x l) inches/millimeters	26 x 36 x 39 660.4 x 914.4 x 990.6	26 x 36 x 48 660.4 x 914.4 x 1219.2	26 x 36 x 60 660.4 x 914.4 x 1524	32 x 36 x 48 812.8 x 914.4 x 1219.2	36 x 42 x 60 914.4 x 1066.8 x 1524	36 x 42 x 84 914.4 x 1066.8 x 2133.6
Chamber Capacity	21 cu. ft. / .59 cu. m.	26 cu. ft. / .74 cu. m.	33 cu. ft. / .93 cu. m.	32 cu. ft. / .91 cu. m.	52.5 cu. ft. / 1.49 cu. m.	73.5 cu. ft. / 2.08 cu. m.
Overall Width	75.00 / 1905	75.00 / 1905	75.00 / 1905	87 / 2209.8	97 / 2463.8	97 / 2463.8
Knockdown Width	37.75 / 958.9	37.75 / 958.9	37.75 / 958.9	43.75 / 1111	49.13 / 1247.9	49.13 / 1247.9
Overall Height ¹	81.00 / 2057.4	81.00 / 2057.4	81.00 / 2057.4	81.00 / 2057.4	86.00 / 2184.4	86.00 / 2184.4
Overall Length (SD) ²	53.19 / 1351.02	62.19 / 1579.62	74.19 / 1884.4	60.5 / 1536.7	73.75 / 1873.3	97.75 / 2483
Wall Opening Width ³	76.00 / 1930.4	76.00 / 1930.4	76.00 / 1930.4	88.00 / 2235.2	100.00 / 2540	100.00 / 2540
Wall Opening Height	80.00 / 2032.0	80.00 / 2032.0	80.00 / 2032.0	80.00 / 2032.0	85.00 / 2159	85.00 / 2159
Wall Opening Length	51.50 / 1308.1	60.50 / 1536.7	72.50 / 1841.5	62.5 / 1587.5	75.75 / 1924	99.75 / 2533.7
Fascia Width ¹	78.00 / 1984.2	78.00 / 1981.2	78.00 / 1981.2	90.00 / 2286	102.00 / 2590.8	102.00 / 2590.8
Fascia Height	81.00 / 2057.4	81.00 / 2057.4	81.00 / 2057.4	81.00 / 2057.4	86.00 / 2184.4	86.00 / 2184.4
Floor to Chamber	24.00 / 609.6	24.00 / 609.6	24.00 / 609.6	24.00 / 609.6	24.00 / 609.6	24.00 / 609.6


1. Fascia extends 1" beyond wall opening on each side and top overlapping the wall and sealing the opening.
2. Allow minimum 2" clearance at rear of recessed unit only. Cabinet side models have rear clearance built in.
3. Contact PRIMUS for alternative wall opening.
4. Standard Left Side Service/equipment access shown. Optional right side access is available.



PRIMUS General Purpose- Medium Sized Autoclaves


SPECIFICATION DATA SHEET

DOUBLE DOOR



SERVICE SIDE

ELEVATION VIEW



PLAN VIEWS

VOLUME/DIMENSIONS CHART

*Refer to General Arrangement (GA) drawings for details and final connection point to utility services (S-Steam, W-Water, D-Drain, E-Electrical, A-Air).

MODEL	E	F	G	G1	H	I
Chamber Size (w x h x l) inches/millimeters	26 x 36 x 39 660.4 x 914.4 x 990.6	26 x 36 x 48 660.4 x 914.4 x 2192	26 x 36 x 60 660.4 x 914.4 x 1524	32 x 36 x 48 812.8 x 914.4 x 1219.2	36 x 42 x 60 914.4 x 1066.8 x 1524	36 x 42 x 84 914.4 x 1066.8 x 2133.6
Chamber Capacity	21 cu. ft. / .59 cu. m.	26 cu. ft. / .74 cu. m.	33 cu. ft. / .93 cu. m.	32 cu. ft. / .91 cu. m.	52.5 cu. ft. / 1.49 cu. m.	73.5 cu. ft. / 2.08 cu. m.
Overall Width	75.00 / 1905	75.00 / 1905	75.00 / 1905	87 / 2209.8	97 / 2463.8	97 / 2463.8
Knockdown Width	37.75 / 958.9	37.75 / 958.9	37.75 / 958.9	43.75 / 1111	49.13 / 1247.9	49.13 / 1247.9
Overall Height ¹	81.00 / 2057.4	81.00 / 2057.4	81.00 / 2057.4	81.00 / 2057.4	86.00 / 2184.4	86.00 / 2184.4
Overall Length (DD)	56.19 / 1427.2	65.19 / 1655.8	77.19 / 1960.6	64.63 / 1641.6	77 / 1955.8	101 / 2565.4
Wall Opening Width ³	76.00 / 1930.4	76.00 / 1930.4	76.00 / 1930.4	88.00 / 2235.2	100.00 / 2540	100.00 / 2540
Wall Opening Height	80.00 / 2032.0	80.00 / 2032.0	80.00 / 2032.0	80.00 / 2032	85.00 / 2159	85.00 / 2159
Wall Opening Length	50.90 / 1292.86	59.90 / 1521.46	71.90 / 1826.26	62.00 / 1574.8	74.19 / 1884.4	98.19 / 2494
Fascia Width ¹	78.00 / 1981.2	78.00 / 1981.2	78.00 / 1981.2	90.00 / 2286	102.00 / 2590.8	102.00 / 2590.8
Fascia Height	81.00 / 2057.4	81.00 / 2057.4	81.00 / 2057.4	81.00 / 2057.4	86.00 / 2184.4	86.00 / 2184.4
Floor to Chamber	24.00 / 609.6	24.00 / 609.6	24.00 / 609.6	24.00 / 609.6	24.00 / 609.6	24.00 / 609.6

1. Fascia extends 1" beyond wall opening on each side and top overlapping the wall and sealing the opening.
2. Contact PRIMUS for alternative wall opening.
3. Standard Left Side Service/equipment access shown. Optional right side access is available.

ARCHITECTURAL NOTES:

1. Allow sufficient space for traps, shut-off's, filters and other utility supply components.
2. Utility (service disconnects) shall be provided and installed "By Others".
3. Provide maximum mechanical and service access space, a minimum of 24", additional space required when boiler specified. See General Arrangement drawing for placement of ancillary equipment and service access.
4. Some options affect utility services and overall dimensions.
5. The Manufacturer's Equipment Warranty does not cover failure due to improper utility provisions.
6. Drawings not to scale.
7. Wall thickness must be provided on double door models recessed through 1 wall, with cabinet sides.

UTILITY SERVICES

Provide utility services within 6'-0" of final connection to sterilizer. Optimum sterilizer performance requires the specified utilities.

MODEL	STEAM (S) Building Steam Supply . Pipe Size: See Below . Quality: Condensate free 97% to 100% saturated vapor (suitably trapped to ensure dry steam and filtered to remove particulates) . Pressure: 50-80 PSIG Dynamic Note: 1.) Steam-to-steam generator requires minimum pressure 65 PSIG house steam		WATER (W) Cold Water Supply . Pipe Size: See Below . Temperature: < 70° F . Pressure: 50-70 PSIG Dynamic		DRAIN (D) Building Drain System Minimum 2" . Location: Locate floor sink in the service area. Note: 1.) Exhaust discharge is cooled to < 140°F 2.) 12" x 12" x 8" floor sink is recommended by PRIMUS 3.) Floor should be sloped to drain	Air (A) Instrument Air . Connection: See Below . Quality: Dry and oil free . Pressure: 60-80 PSI Dynamic	ELECTRICAL (E) Building Power Supply . Volts: 110 . Phase: Single . Amps: 10 Note: Provide dedicated circuit to sterilizer controls. If optional boiler is specified, do not share this service
	NPT	LBS/HR (KG/HR)	NPT	GPM (MAX) (Liter)	NPT (Discharge Pipe Size)	NPT	
E	3/4"	169.0(76.66)	3/4"	14 (53)	1"	1/4"	
F	3/4"	200.0 (90.72)	3/4"	14 (53)	1"	1/4"	
G	1"	260.0 (117.94)	3/4"	14 (53)	1"	1/4"	
G.1	1"	256.0 (116.12)	3/4"	14 (53)	1"	1/4"	
H	1"	350.0 (58.76)	1"	14 (53)	1 1/4"	1/4"	
I	1"	420.0 (190.51)	1"	25 (76)	1 1/4"	1/4"	

HVAC DATA Heat loss, at ambient of 70° F.

	MODEL	KBTU'S/HR
SINGLE DOOR: Through one wall, at fascia	E	5
	F	5
	G	5
	G.1	5.2
	H	6.8
	I	6.8
SINGLE DOOR: Through one wall, service area	E	7.5
	F	9
	G	11.8
	G.1	12.2
	H	10.5
	I	12.3
SINGLE DOOR: Free standing, cabinet total	E	13.3
	F	14
	G	16
	G.1	16
	H	16.1
	I	16.3

	MODEL	KBTU'S/HR
DOUBLE DOOR: Through one wall, at fascia	E	5
	F	5
	G	5
	G.1	5.2
	H	6.8
	I	6.8
DOUBLE DOOR: Through one wall, service area	E	12.3
	F	13
	G	16
	G.1	13.8
	H	15
	I	16.3
DOUBLE DOOR: Through two walls, at each fascia	E	5
	F	5
	G	5
	G.1	5.2
	H	6.8
	I	6.8
DOUBLE DOOR: Through two walls, service area	E	6.9
	F	8
	G	10
	G.1	9.2
	H	7.9
	I	10.2

Use the following pages to custom design your sterilizer by checking the boxes of the configuration and options required for your project. Transfer all items checked to the Specification Worksheet (Pages 11 & 12) and fax to PRIMUS Sales Department or to your local Sales Representative for a quotation. If you do not find an option or size listed to meet your requirements, please contact our Project Manager for additional information and assistance.

Model Number Notation

SD = Single Door
DD = Double Door
↓

PSS - -M - -

↑ ↑ ↑

PSS5 Control = 5	E	S = House Steam	Five models are offered:
Pri-Matic = 6	F	E = Electric Boiler	E - 26" x 36" x 39"
ControlLogix = 7	G	P = Steam to Steam	F - 26" x 36" x 48"
Trinity Control = 8	G.1		G - 26" x 36" x 60"
	H		G.1 - 32" x 36" x 48"
	I		H - 36" x 42" x 60"
			I - 36" x 42" x 84"

A detailed formatted specification, suitable for inclusion in formal contract documents, is available on request.

Note: Use the above format to determine PRIMUS Model number and insert below and on Page 11, Specification Worksheet.

PRIMUS General Purpose Steam Pressure Sterilizer, Model Number PSS_____.

SELECT CONFIGURATION

Door

- GDA - Single Door
- GDB - Double Door

Cabinet

- GCD - Panels Both sides
- GCCL - Left Side Panel
- GCCR - Right Side Panel

Recessed

- GCA - One wall
- GCB - Two walls

Service and Equipment Access

- Left Side (Standard)
- Right Side

Steam Source

- | | | |
|--|---|---|
| <input type="checkbox"/> House Supply | <input type="checkbox"/> Clean Steam ² | <input type="checkbox"/> Steam to Steam ² |
| <input type="checkbox"/> Electric Boiler ^{1,2,3} | <input type="checkbox"/> GEBC1 - 208/3 phase ³ | <input type="checkbox"/> GSCS1 ⁴ - Pure Steam Dbl Tube |
| <input type="checkbox"/> GEB1 - 208/3 phase | <input type="checkbox"/> GEBC2 - 240/3 phase | <input type="checkbox"/> GSCS2 - Clean Steam Sgl Tube |
| <input type="checkbox"/> GEB2 - 240/3 phase | <input type="checkbox"/> GEBC3 - 480/3 phase | <input type="checkbox"/> GSCS3 ⁴ - Clean Steam Sgl Tube Sanitary |
| <input type="checkbox"/> GEB3 - 480/3 phase | <input type="checkbox"/> Other _____ | <input type="checkbox"/> Other _____ |
| <input type="checkbox"/> Other _____ | | |
| <input type="checkbox"/> GEB0 ⁵ - Auto Blowdown | | |

1. Electric boilers are available in Carbon Steel or Stainless Steel.
2. Contact PRIMUS for overall dimensions and utility connections.
3. All models are stand-alone.
4. GMP Validatable Option Only, sanitary piping required.
5. If Option LWC is required, "the automatic" feature is disabled. The boiler will need to be manually reset.

RESEARCH & LABORATORY CONFIGURATIONS

Clean Steam Sterilizers GP9

For cell/tissue culture applications. This option provides direct connection to a house-source of clean steam. Includes threaded stainless piping and components to chamber. If clean steam is the only steam source available, specify option GP9.1 for stainless piping to jacket. For alternate sources to generate clean steam, specify electric clean steam boilers (stainless steel), or when house steam is available, specify steam-to-steam generators. (Stainless steel heat exchanger).

GP9	Clean Steam to Chamber, SS Threaded Piping
GP9.1	Clean Steam to Jacket, SS Threaded Piping

Note: See Steam Source on previous page

Low Temperature Flowing Steam GP14

For media preparation, formula and similar uses. Pasteurization is accomplished with flowing steam. In addition to standard sterilizer temperature ranges,

this option allows cycles to operate in the range of 168°F (76°C) to 220°F (104°C). Not available when effluent configured sterilizer option GP8 is specified.

Note: GP14 not available if option GP8, Effluent Configured Sterilizer, is specified.

Air-over Cooling GP11

Provides a means of rapidly cooling sterilized load. Effective to decrease cycle times with liquid loads. Filtered ambient temperature air displaces steam under pressure during exhaust cycle to rapidly cool load and prevent boil-over.

Control Cycle on Load Temperature GC24

Used to establish the most beneficial cycle temperature and exhaust rate parameters. Monitors load temperature and controls the start of sterilization. Which will establish optimal parameters for liquid loads and for loads with complex load geometries and a difficult path for steam penetration.

Note: See additional Bio-Containment Options

GMP VALIDATABLE CONFIGURATIONS

PRIMUS offers our Bio→Pharma® line of validatable sterilizers for Pharmaceutical applications. PRIMUS General Purpose sterilizers may be configured to be validatable, in accordance with the customer's specific requirements. To assure a compatible configuration is being specified, complete the worksheet and submit to PRIMUS Validatable Sterilizer Product Manager for review at validatable@primus-sterilizer.com.

CONTROLS AND RELATED OPTIONS

(Contact PRIMUS for details and performance specifications)

PSS500 Control GC6

PSS500 microcomputer with eight (8) cycles; seven (7) programmable and 1 preset Test cycle are standard. A vertical touchpad contains a Liquid Crystal Display (LCD) mounted on a side-mounted control panel with a thermal printer. Cycle progression is displayed by Light Emitting Diodes (LEDs) on the operator panel and printed, 32 characters per line printing. A second copy of the entire sterilization cycle printable at the end of the run. (If the reprint button is pressed on the printer, Cybertech printer only.)

PSS6 Pri-Matic® Control GC1

Pri-Matic® Control System uses the Allen-Bradley SLC 500 Programmable Logic Controller (PLC) with an Allen-Bradley Panel View 1000 Operator Interface Terminal (OIT) and the Pri-Matic® software package to provide automatic cycle operation. Pri-Matic® offers a total of fifteen cycles, with thirteen standard cycles available for product sterilization. Standard cycles are configurable by modifying the cycle parameters for each. Two dedicated cycles are available for performing vacuum leak test and sterilize filter.

PSS7 Control System GC1.1

PSS7 uses the Allen-Bradley ControlLogix Series 5555 PLC (Ethernet) and PanelView Plus OIT in place of the SLC500 PLC and standard PanelView. ControlLogix provides graphical color touch screen capability, real time process graph and the ability to standardize control systems throughout a facility. PSS7 offers a total of fifteen cycles, with thirteen standard cycles available for product sterilization. Standard cycles are configurable by modifying the cycle parameters for each. Two dedicated cycles are available for performing vacuum leak test and sterilizer filter.

PSS8 Trinity Control GCT

PSS800 microcomputer with a selection of twenty-five(25) separate programs for sterilization of wrapped goods, hard-goods, liquids and test cycles. Trinity controller has a color touch sensitive screen with 30 line x 40 character display on a side-mounted control panel with a thermal printer. Help screens for programming and trouble shooting alarm conditions. Cycle data can be printed, captured to a remote personal computer or logged to memory.

CONTROLS AND RELATED OPTIONS (Con't)

(Contact PRIMUS for details and performance specifications)

Serial Data Output GC11
Supervisory use to capture electronic record of each sterilization cycle. Cycle hard copy and display data is transmitted via RS232 in real time to a remote computer, up to 100 feet away, for display and recording. Hard copy records may be printed from the computer. Available on both control systems.

Authorized User Access GC32
Restricts operation of sterilizer to authorized personnel and provides supervisor with PIN number of person running the cycle.

Remote Mount Control Panel GC10
Mount control panel in separate housing adjacent to or up to 100 feet away from the sterilizer. Available on both control systems.

Display & Record in PSIA GC25
Pressure values are displayed and recorded in pounds per square inch absolute in lieu of standard gauge pressure (PSIG).

Thermal Printer GR7
Thermal dot-matrix printer with take up reel and 32 character per line printing is standard. Second sterilization cycle report available at the end of the run.

BIO-CONTAINMENT OPTIONS

Effluent Controlled Configuration GP8
For use in high-risk biohazard applications including BSL-3 laboratories. All effluent is either filtered or contained within the chamber until completion of the sterilization phase of the cycle. In addition to a specially designed chamber, GP8 option includes special components making it suitable for biohazard application. *Contact PRIMUS for details and performance specifications.*

BioSeal Flange GV6
Provides a means of isolating load from the unload ends of the sterilizer. A vapor-proof flange is welded to the vessel on both sides of the flange and around its full circumference. Penetrations through the flange for electrical and piping components are secured with vapor proof fittings.

BioSeal Panels GV7
Aluminum panels for attachment to, and extending from, the BioSeal Flange to adjacent building surfaces to complete the barrier between either end of the sterilizer.

	GV7	BioSeal Extension Panels, Aluminum
	GV10	BioSeal Extension Panels, SS

Compressed Air to Gasket GP10
Compressed air, in lieu of steam, is used to seal gasket against door interior during cycles. Prolongs gasket life. Requires compressed air utility service or option GEAC air-compressor.

Remote Signaling of Sterilizer Status GC23
Control relay connection enabling operation of a remote signaling devise for an alarm condition. Alarm conditions many include Door Open, Door Closed and Cycle complete. Signal may be used to activate a buzzer, light or any other On-Off device.

GENERAL OPTIONS

Validation Port, 1.5 inch GV4
Provides chamber penetration to accommodate various monitoring/control probes. Standard on all models.

Seismic Restraints GV8
Required in areas prone to seismic hazards. Secures sterilizer to building but allows for leveling. Designed to current California Code.

Electric Vacuum Pump GP4
This water-sealed electric vacuum pump provides an alternative means of drawing a vacuum, in lieu of the standard water ejector. Check PRIMUS' General Arrangement drawing to assure the required ancillary equipment space is provided.

Drain Line Strainer/Valve GP25
Provides a means to trap and expel debris from the chamber drain line to protect heat exchanger from damage.

Air Compressor GEAC
For sterilizers, less than 75 cu. ft., air compressor mounted on sterilizer frame, for use when house supplied air is required but not available. Not available if Air-Over Option GP11 is specified.

UTILITY - RELATED OPTIONS

Boost Pump for Low House Water Pressure GP3

Delivers required dynamic water pressure for efficient operation. Requires specified GPM flow. *Contact PRIMUS for details and performance specifications.*

Step Down Transformers GE1

Reduces line voltage to required 110 VAC for operation of sterilizer controls.

	GE1A	220 VAC to 110 VAC
	GE1B	480 VAC to 110 VAC
	GE1C	480/240 to 240/120 VAC, 1 phase, 0.5 KVA

Uninterruptible Power Supply GE2

In the event of electrical power loss and no emergency electrical power to sterilizer, this option provides electrical power to operate control system for up to 30 minutes, to complete the cycle.

Water Recirculation Conservation System GP24

Provides water conservation through the use of Chilled Water versus City Water to cool drain effluent. Additional water conservation is provided by reusing the cooled drain effluent for the water ring seal on the vacuum pump.

LOADING EQUIPMENT

Standard chamber shelving includes a fixed, removable, wire mesh bottom shelf. As an option an extendable bottom shelf, is available. Optional additional chamber shelf(ves) may be specified.

Loading cart includes one bottom and one intermediate shelf with four adjustable levels. Additional shelves are available. Cart frames and shelves are 316L stainless steel welded, ground and polished. Shelf surfaces are stainless steel wire mesh. Transfer carriages include swivel casters with swivel locks and 5" wheels with wheel brakes. Loading cart and transfer carriages are recommended for chamber lengths 48 inches or longer.

	GL1	Removable Bottom Shelf
	GL2	Extendable Bottom Shelf (26"x36"x39" Only)
	GL3	Additional Chamber Shelf (26"x36"x39" Only)
	GL4	Loading Cart
	GL5	Additional Cart Shelf
	GL6	Transfer Carriage

Contact PRIMUS for special loading equipment requirements.

SHIPPING DIMENSIONS, CUBAGE & WEIGHTS

Model Sterilizer Size Inches/millimeters	E 26" x 36" x 39" 660.4 x 914.4 x 990.6	F 26" x 36" x 48" 660.4 x 914.4 x 1219.2	G 26" x 36" x 60" 660.4 x 914.4 x 1524	G1 32" x 36" x 48" 812.8 x 914.4x 1219.2	H ¹ 36" x 42" x 60" 914.4 x 1066.8 x 1524	I ¹ 36" x 42" x 84" 914.4 x 1066.8 x 2133.6
Shipping Dimensions, Single Door ²	75 x 70.25 x 48.5 1905 x 1784.4 x 1231.9	75 x 70.25 x 57.5 1905 x 1784.4 x 1460.5	75 x 70.25 x 69.5 1905 x 1784.4 x 1765.3	87 x 70.25 x 58.25 2209.8 x 1784.4 x 1479.6	65.5 x 77.88 x 71.5 1663.7 x 1978 x 1816.1	65.5 x 77.88 x 95.5 1663.7 x 1978 x 2426
Shipping Dimensions, Double Door ²	75 x 70.25 x 49 1905 x 1784.4 x 1244.6	75 x 70.25 x 58 1905 x 1787.4 x 1473.2	75 x 70.25 x 70 1905 x 1784.4 x 1778	87 x 70.25 x 60 2209.8 x 1784.4 x 1524	65.5 x 77.88 x 72.38 1663.7 x 1978 x 1838.3	65.5 x 77.88 x 96.38 1663.7 x 1978 x 2448
Overall Width Dimensions	75 / 1905	75 / 1905	75 / 1905	87 / 2209.8	97 / 2463.8	97 / 2463.8
Knockdown Width Dimensions ³	37.75 / 958.9	37.75 / 958.9	37.75 / 958.9	43.75 / 1111.3	49.13 / 1247.9	49.13 / 1247.9
Weight, Single Door (with plumbing)	CONTACT PRIMUS FOR WEIGHTS					
Weight, Double Door (with Plumbing)						
Crated Weight (Additional)						
	VAN DELIVERY ONLY					
Crated Dimensions	78 x 80 x 58 1981.2 x 2032 x 1473.2	78 x 80 x 67 1981.2 x 2032 x 1701.8	78 x 80 x 79 1981.2 x 2032 x 2006.6	90 x 80 x 69 2057.4 x 2032 x 1752.6	81 x 83 x 86 2057.4 x 2108.2 x 2184.4	81 x 83 x 110 2057.4 x 2108.2 x 2794
Crated Cube	209 cu ft / 5.9 cu m	242 cu ft / 6.9 cu m	285cu ft/ 4 cu m	259 cu ft/ 7.3 cu m	335 cu ft/ 9.5 cu m	428 cu ft/ 12.1 cu m
	OPTIONAL ACCESSORIES					
Loading Equipment	146 lbs / 66.kg	152 lbs/ 69kg	205 lbs/ 92.99 kg	166 lbs/ 75.3 kg	282 lbs/ 128 kg	370 lbs/ 168 kg
Boiler	380 lbs/ 172 kg	390 lbs/ 177kg	625 lbs/ 283.5 kg	625 lbs/ 284 kg	785 lbs/ 356 kg	785 lbs/ 356 kg

- Units have Split Beams and may be split for shipping purpose. All other units do not have split beam.
- Shipping dimensions are measured to the edge of the heat exchanger and include plumbing on, split beams removed, where applicable, front stainless steel panel and printer removed.
- Knockdown is vessel dimension only, dimension does not include plumbing, electrical or heat exchanger. Special Order Only, must be indicated when ordering.

***“PRIMUS Sterilizer...
...First In Quality Sterilization Products”***



***26"x36"x48" with loading
cart
and transfer carriage***

***QUALITY AND INTEGRITY FORM THE FOUNDATION OF
PRIMUS STERILIZER CO., INC.
OUR VISION FOR CONTINUAL IMPROVEMENT IS REALIZED
THROUGH THE ESTABLISHMENT OF QUALITY OBJECTIVES,
ACHIEVED THROUGH TEAMWORK WITH A GOAL OF
RETURNING SIGNIFICANT VALUE TO OUR CUSTOMERS.***



Specification Data Worksheet (con't)

LOADING EQUIPMENT	
GL1	Removable Bottom Shelf (Std)
GL2	Extendable Bottom Shelf (26"x36"x39" Only)
GL3	Additional Chamber Shelves (26"x36"x39" Only)
GL4	Loading Cart
GL5	Additional Cart Shelves
GL6	Transfer Carriage
OTHER CONTROL OPTIONS	
GC2 ³	Pri-Matic® F ₀ control F ₀ accumulated, 4 points
GC5 ³	Conax Adaptor
GC7 ^{2,3}	Second PanelView 1000 w/ DH+
GC8 ^{2,3}	SCADA Data File
GC9 ¹	Display Centigrade
GC12	Chamber Condensate Alarm
GC13 ^{2,3}	Sterilize Filter as a Cycle.
GC14 ^{2,3}	WFI Interface
GC15 ^{2,3}	Pri-Matic® in NEMA 12 enclosure (Dustproof)
GC15.1 ³	Stainless Steel Enclosure for Pri-Matic Control Panel
GC16 ^{2,3}	Pri-Matic® mounted in NEMA 4 Enclosure (Waterproof)
GC16.1 ^{2,3}	Pri-Matic® mounted in NEMA 4X SS Enclosure
GR1	Impact Printer, Flush Mount 24 Col.
GR3 ³	Strip-Chart Recorder
GR4	Thermal Printer, Take-up Reel
GR4.1	Thermal Printer, Flush Mount 24 Col.
GR5 ³	Circular Chart Recorder
GR6	Impact Printer, 24 Column

- 1. PSS500 Control Only
- 2. Pri-Matic® Control Only
- 3. Validatable Option Only
- Std = Standard Component
- SS = Stainless Steel
- CIP = Clean-in-Place
- N/A = Not Available

OTHER CONTROL OPTIONS (con't)	
GC17 ³	Sterilizer "OFF" Signal
GC18 ³	Sterilizer "In-Cycle" Signal
GC19	Chamber Drain Temperature, Remote Digital Display
GC23 ^{2,3}	Remote Alarm Signal
GC26 ^{2,3}	PanelView 550, Service Area
GC27 ^{2,3}	PanelView(s) Color
GC28 ^{2,3}	Ethernet PLC
GC33	Automatic Vacuum Leak Rate Test

OTHER PIPING OPTIONS	
GP1	Heated Air-In, Heat Exchanger SS
GP2	Heated Air-In, Heat Exchanger Non-SS
GP5 ³	Safety Valve Over Rupture Disk
GP7 ³	Condensate Sampling Valve
GP13	Precise Temperature Control +/- .5°C (S3)
GP15 ³	CIP Piping for Sanitary Chamber
GP16 ^{2,3}	Air to Jacket Cooling
GP18 ³	Passivate Chamber
GP19 ³	Passivate Chamber and Piping
GP20 ³	Steam Quality Sampler
GP21 ³	Adapters, Sanitary to Threaded
GP24	Water Recirculation Conservation System

VESSEL OPTIONS	
GV1	Stainless Steel Frame
GV2	Jacket Type, 316 L SS
GV3	Jacket Insulation Cover, SS
GV6.1	Additional BioSeal Flange
GV8	Seismic Restraints
GV9 ³	Chamber Welds Ground and Polished
GV10	BioSeal Panels, Stainless
GV12 ³	Additional Validation Port
GC3 ³	Blind Flange Assembly, Sanitary fittings
GC4	Jacket Pressure Display, Analog

GMP DOCUMENTATION OPTIONS	
GD1 ³	GMP Validation Documentation
GD2 ³	IQ/OQ Protocol, PSS500
GD3 ³	IQ/OQ Protocol, Pri-Matic®

YOUR LOCAL PRIMUS REPRESENTATIVE:

FOR ADDITIONAL ASSISTANCE

To discuss your sterilization design needs or to request a Sterilizer Design Guide, please contact TES.

Manufacturer reserves the right to modify materials or specification without notice.

